TheraMart

TM2501

TM2520

PT

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Cadeira de Rodas Motorizada

ES

MANUAL DE INSTRUCCIONES

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EN

INSTRUCTION MANUAL

Motorized Wheelchair

Prefácio

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[EN]

Preface

Please read the user manual carefully before using the product. This manual provides instructions on operation, assembly methods, and basic troubleshooting for the following models: **TM2501, TM2520.**

This manual includes wheelchair maintenance tips and troubleshooting guidelines. Please keep it in an accessible location and refer to it as needed. If someone else uses the wheelchair, please provide them with this manual for reference.

The notes and illustrations in this manual may differ slightly from your actual model, as we continuously improve and update our components and designs.

If you are unsure about any aspect of the operation, please contact your dealer for clarification.

Improper use of this wheelchair may result in injury. It is essential to always follow safety recommendations to ensure your well-being and that of others.

This electric wheelchair is designed for individuals with mobility difficulties or those who cannot walk.

The electric wheelchair is intended to carry only one person.

This series of electric wheelchairs is powered by lithium-ion batteries and uses DC motors. Users can control the direction and adjust the speed using the joystick controls located on the armrest. This wheelchair is suitable for low speeds, good road conditions, and can handle gentle inclines.

This electric wheelchair is specifically designed for people weighing less than the maximum load weight set in the technical specifications. It is intended for walking speeds and is suitable for both indoor and outdoor use. Please avoid using this unit in traffic lanes.

This user manual is compiled based on the latest specifications and product information available at the time of publication. We reserve the right to make changes as necessary. Any changes to our products may result in slight variations between the illustrations and explanations in this manual and the product you purchased.

Revision 03

Date: December 2024

Warning Symbol

Carefully follow the instructions next to this symbol.

Please pay attention to these instructions; otherwise, it may result in physical injury or damage to the wheelchair or the environment.

1. Backrest
2. Controller (Joystick)
3. Controller Charging Port
4. Safety Belt
5. Cushion (Seat Base)
6. Shopping Bag
7. Foot Platform (Footrest)
8. Front Wheel (Support Wheel)
9. Rear Wheel (Drive Wheel)
10. Lithium Battery

11. Armrest Bag
12. Armrest Cushion
13. Backrest Cushion
14. Backrest Cushion
15. Rear Bag
16. Folding and Locking Mechanism (Safety Lock)
17. Motor
18. Manual Unlock Lever
19. Anti-Tip Wheel

Product Performance

This series of electric wheelchairs operates with lithium-ion batteries and is powered by a direct current (DC) motor.

Users control the direction and adjust the speed using the joystick control.

It is suitable for low-speed applications, good road conditions, and small inclines.

Product Information

Performance Parameters

Main Technical Data	TM2501	TM2520
		
Total Size	98x60x95 cm	98x60x95 cm
Seat Size	45*44 cm	45*44 cm
Speed	1~6 km/h	1~6 km/h
Armrest Height	24 cm	24 cm
Maximum Load Capacity	100 kg	100 kg
Rear Wheel Diameter	12 inches	12 inches
Motor	DC 24V 250W*2	DC 24V 250W*2
Control	24V 50A	24V 50A
Obstacle Clearance	5 cm	5 cm
Turning Radius	<1.2 m	<1.2 m
Turning Width	<1.5 m	<1.5 m
Tire Pressure	2kgf/m ²	2kgf/m ²
Folded Size	80x68x40 cm	80x68x40 cm
Net Weight (N.W.)	27/28 kg	27/28 kg
Seat Height	52 cm	52 cm
Backrest Height	50 cm	50 cm
Front Wheel Diameter	8 inches	8 inches
Battery	DC 24V 12Ah/20Ah Li-ion	DC 24V 12Ah/20Ah Li-ion
Battery Charger	AC 100~240V, 50Hz, 1.5~5A	AC 100~240V, 50Hz, 1.5~5A
Maximum Driving Distance	15/25 km	15/25 km
Maximum Safe Inclination	10°	10°
Static Stability	9°	9°
Climbing Angle	6°	6°

Headrest	No	Yes
Reclining Backrest	No	Yes

Application Area

Our electric wheelchair is designed for people with disabilities and elderly individuals who use it as a means of mobility, suitable for indoor use. It should not be used on roads or public streets.

Contraindications

The following individuals should not use this wheelchair: Lack of response in upper limbs, Senile dementia, Psychiatric disorders, Physiological inability for self-care, Medical recommendation against its use

2. Safety Instructions

2.11 Main Safety Features

- Classification according to protection type against electric shock: Internal power supply
- Classification according to safety in an inflammable anesthetic mixture with air or gas mixed with oxygen or nitrous oxide: Not AP/APG type
- Classification according to operation mode: Continuous operation
- Rated voltage: DC 24V
- Does not provide protective effect against defibrillator discharges
- No signal input or output
- Equipment is not permanently installed

2.12 Driving Warnings

General Driving Precautions

- Keep your hands on the armrests to control the wheelchair.
- Practice in parks or open and safe areas until you can handle the wheelchair skillfully.
- Familiarize yourself with driving in a safe environment before using it in other spaces.
- Before driving on public roads, ensure you are accompanied by a caregiver and check that it is safe.
- Strictly follow pedestrian traffic rules. Do not assume you are driving a vehicle.
- Drive only on sidewalks and pedestrian crossings. Do not drive in vehicle lanes.
- Avoid sudden or zigzag movements that may compromise stability.
- Maintain proper tire pressure to prevent instability and excessive battery consumption.

Conditions in which you should be accompanied by a caregiver or avoid driving

- Driving in adverse weather conditions such as rainy days, dense fog, strong winds, snow, etc.
- If the wheelchair gets wet, dry it immediately.
- Driving on poorly maintained roads such as mud, trails, sand paths, gravel, etc.
- Driving on congested streets.
- Driving near unprotected ditches, lakes, etc.
- Crossing railway tracks.
- When crossing a railway track, stop at the crossing to ensure it is safe and check that the tires do not get stuck on the rails.
- The electric wheelchair is for personal use only. Do not transport people or goods and do not use it for towing.

Precautions for Going Up and Down Slopes

- Avoid driving in the following places: steep hills, inclined surfaces, high steps, channels, etc.
- Avoid driving on steep slopes; the incline should be less than 9 degrees.
- Operate the controller carefully when driving on slopes.
- Continue moving forward while ascending or descending.
- Reduce speed when descending.
- Avoid driving sideways on steep slopes.
- It is prohibited to drive over steps and avoid crossing high steps.
- Avoid crossing wide ditches.
- When crossing ditches, maintain a 90° angle between the tires and the ditch.
- Do not place the wheelchair in manual mode while going up or down slopes.

If the wheelchair malfunctions at a traffic crossing, immediately request assistance from passersby. Then, set the wheelchair to manual mode, push it out of the area, or move the user to a safe location.

Precautions for Caregivers

Caregivers must ensure the user's feet are correctly positioned on the footrests and make sure that clothing does not get caught in the wheels.

Caregivers should push the wheelchair forward to maintain safety on steep or prolonged slopes.

2.13 Other Considerations

Repair and Modification

If the wheelchair requires repair or modification, contact the seller or service department. Do not perform modifications on your own, as this may cause accidents or malfunctions in the wheelchair.

Keep Dry

Do not place the wheelchair in humid areas. If the wheelchair gets wet, dry it immediately.

Without our company's confirmation, do not modify the assembly or materials of the wheelchair. To avoid imbalance, do not arbitrarily add weight. When someone is seated in the wheelchair or when

the clutch is not in manual mode, do not use other vehicles to push or tow the wheelchair.

2.14 Lifespan

The lifespan of this product is 5 years from the date of manufacture. Please use the product within the specified period. Do not use it beyond the validity period for more than one year to avoid accidents. Date of Manufacture: Please refer to the product label.

2.15 Electromagnetic Compatibility

The wheelchair must be kept away from strong magnetic fields and high-power electrical equipment, such as: Radio stations, Television stations, Underground radio stations, Cell phone signal transmission towers.

If you detect sources of electromagnetic interference nearby, move away from them to avoid disruptions in the wheelchair's operation. The electric wheelchair must avoid electromagnetic interference.

Note:

- The electric wheelchair must comply with the electromagnetic compatibility requirements of standard YY0505.
- The user must install and use the electric wheelchair according to the electromagnetic compatibility information provided.
- Portable and mobile RF communication devices may affect the performance of the electric wheelchair.
- Therefore, avoid strong sources of electromagnetic interference, such as mobile phones or microwaves.
- Refer to the annex for the manufacturer's warning and declaration.
- The battery cut-off voltage is 23V.
- The electric wheelchair belongs to Class D according to the GB/T 18029.21-2012 standard, with differential electronic steering and manual braking.

The electric wheelchair should not be placed or used together with other devices. If necessary, verify that the electric wheelchair can operate normally under those conditions.

3 Use and Operation

3.1 Wheelchair Adjustment

Unfolding the Wheelchair

Using one hand, hold the backrest and with the other, grasp the seat and open it firmly (Image 1). Fully unfold the wheelchair and ensure that the lock is securely engaged (Image 2).

Picture 1

Picture 2

Ensure the Lock is Fully Engaged. If the lock is not completely secured, there is a risk of accidental folding while driving.

Folding the Wheelchair

Loosen the lock (Image 3), then, with one hand, hold the backrest, and with the other, pull the seat and fold it (Image 4).

Picture 3

Picture 4

Installing the Controller

Insert the controller into the armrest tube. Make sure to keep the controller in a horizontal position (Image 5). Then, tighten the screw to secure it (Image 6).

Picture 5

Picture 6

Operation of the Versions with Automatic Reclining and Automatic Folding (TM2520)

Automatic Reclining: If your wheelchair has an automatic reclining function, note that it cannot be reclined manually.

Removing the Battery

Disconnect the controller plug from the battery connector (Image 7).

Loosen the screw and then remove the battery (Image 8).


Picture 7

Picture 8

Installing the Battery

Place the battery groove into the rear tube. Install the screws and nuts (Image 9). Then, connect the controller plug to the battery connector (Image 10).

Picture 9

Picture 10

Footrest Platform

Adjusting the footrest extension: The footrest offers three extension options (Image 11).

Adjusting the footrest height: The footrest offers three height options (Image 12).

Picture 11

Picture 12

Ensure the footrest platform is properly secured to avoid safety issues.
The adjustable foot pedal is an optional feature and may not be available from all distributors.

Installing and Removing the Anti-Tip Device

Press the button to install or remove the anti-tip device. The anti-tip system offers two position length options (Images 13-14).

Picture 13

Picture 14

Set it to the longest position when facing a slope.

Flip-Up Armrest

With one hand, press the armrest button, and with the other, lift the armrest (Image 15). Then, press the armrest down to the appropriate position to secure it (Image 16).

Picture 15

Picture 16

Switching Between Manual Mode and Electric Mode

There are small red levers on the left and right rear wheel motors. Pull the levers down to activate manual mode (unlocked) (Image 17). Pull the levers up to activate electric mode/driving (locked) (Image 18). For proper operation, both levers must always be in the same position.

Picture 17

Picture 18

Safety Belt (Positioning Belt)

The wheelchair may be equipped with a safety belt, which can be adjusted for user comfort. The safety belt is designed to prevent the user from sliding forward or down in the seat. It is not designed as a restraint device.

Picture 19

4 Battery and Charger

4.7 Charger Requirements

The charger is used to charge the battery. Do not use the electric wheelchair while it is charging.

Technical Specifications of the Charger:

Input Voltage: AC 220V (110V) $\pm 10\%$

Output Voltage: 24V/DC

Output Current: 1.5~5A

Protection Level: IPX1

The charger must comply with the requirements of GB 4706.1-2005 and GB 4706.18-2005 standards.

4.8 Using the Charger

4.8.1 Type 1 (Direct Charging to the Battery)

To charge the battery, connect the charger to the power supply and to the battery box socket. Follow these instructions to complete the charging process:

1. Ensure that the charger's ventilation slots are not blocked.
2. Verify that the electric wheelchair is turned off.
3. Disconnect the plugs connecting the battery box and the controller.
4. Connect the charger's output plug to the battery socket.
5. Plug the charger's main plug into the power supply.
 - a. The red light will turn on.
 - b. Full charging takes between 8 to 10 hours. Do not charge for more than 24 hours.
6. To avoid reducing battery lifespan, charge it at least once a month if the wheelchair is not in use.

4.8.2 Type 2 (Charging Through the Controller)

To charge the battery through the controller, follow these steps:

1. Ensure that the charger's ventilation slots are not blocked.
2. Verify that the electric wheelchair is turned off.
3. Ensure that the battery and controller are connected.
4. Connect the charger's output plug to the socket underneath the controller.
5. Plug the charger's main plug into the power supply.
 - a. The red light will turn on.
 - b. Full charging takes between 8 to 10 hours. Do not charge for more than 24 hours.
6. To avoid reducing battery lifespan, charge it at least once a month if the wheelchair is not in use.

Do not interrupt the charge until the process is completed. Using the battery without fully charging it will shorten its lifespan, so always charge the battery completely whenever possible. When the battery is fully charged, the power indicator will turn green. Do not stop charging before it is fully completed. After finishing the charge, disconnect the power supply, or the battery will slowly discharge. Do not charge for more than 24 hours.

Overcharging is Dangerous. Users must follow these rules to avoid hazards during charging:

The electric wheelchair does not include a charger; therefore, use a national standard charger with an output voltage of 24V/DC 1.5A~5A.

Ensure proper ventilation during charging. Do not expose the wheelchair to direct sunlight or humid environments.

Recommended ambient temperature for charging: 10°C to 50°C. If outside this range, the battery may not function properly and could be damaged.

It is normal for the fan to make noise during charging, as its function is to cool the charger. Do not worry about this sound.

Avoid allowing liquids to enter the charger during charging. Do not place the charger on flammable items, such as fuel, footrests, or seat cushions.

Keep away from flames while charging the battery, as they can cause fire or explosion. Charging generates hydrogen, so smoking is prohibited during the charging process.

Do not unplug the power source if the plug or your hands are wet, as this may cause an electric shock. In case of an unexpected accident that may endanger the user, do not use or sit on the electric wheelchair while it is charging.

4.9 Battery Usage and Maintenance

Incorrect battery replacement may cause an explosion hazard. Only use the same battery type or the manufacturer's recommended model. Ensure that the battery terminals are properly connected.

Key points to extend battery lifespan: Charge the battery frequently to keep it fully charged. If the wheelchair is not in use, it is best to fully charge the battery. If the wheelchair is unused for a long period, charge it at least twice a month.

5 System Diagnostics

The LED indicators on your device provide signals regarding various technical malfunctions, which may affect components such as: Motor, Brakes, Battery, Wiring Connections. The following table will help you troubleshoot specific technical issues.

Although these problems can usually be identified without additional tools, if you need further assistance, do not hesitate to contact us.

Refer to Table 1-1 for interpretations of different audio signals.

Sound Count	Interpretation	Problem/Solution
One short beep	The battery power is very low.	Please, charge the battery.
One short beep, the red battery indicator is flashing	Please, charge the battery.	Please, charge the battery.
Two short beeps	Malfunction in the left motor circuit.	1. The left motor connection may be broken, or the cable may be disconnected. 2. Please swap the connection harnesses of the left and right motors; if there are four alarms, it is a defect in the motor circuit. If there are two alarms, it is a defect in the controller circuit.

Three short beeps	The left magnetic brake is in manual mode.	1. The left magnetic brake is in manual mode; please switch to electric mode. 2. Please check if the magnetic brake connection harness is intact. 3. Swap the connection harnesses of the left and right motors; if there are five alarms, it is a defect in the left motor brake. If there are three alarms, it is a defect in the brake control circuit.
Four short beeps	Malfunction in the right motor circuit.	1. The right motor connection may be broken, or the cable may be disconnected. 2. Please swap the connection harnesses of the left and right motors; if there are four alarms, it is a defect in the motor circuit. If there are two alarms, it is a defect in the controller circuit.
Five short beeps	The right magnetic brake is in manual mode.	1. The right magnetic brake is in manual mode; please switch to electric mode. 2. Please check if the magnetic brake connection harness is intact. 3. Swap the connection harnesses of the left and right motors; if there are three alarms, it is a defect in the left motor brake. If there are five alarms, it is a defect in the brake control circuit.
Six short beeps	Motor overload.	1. Motor overcurrent protection alarm. 2. When the alarm sounds, the gear indicator light will show 1st gear, indicating overcurrent in the left motor, and 5th gear, indicating overcurrent in the right motor. 3. Please disconnect the left and right motor cables from the controller, restart the controller, and immediately push the joystick. If there is no alarm or if alarms 2, 3, 4, or 5 appear, it indicates a short circuit in the motor. If there are six alarms, it indicates a short circuit in the controller circuit.
Seven short beeps	The joystick is not in the initial position.	1. Turn off the power, release the joystick to the central position, and turn the power back on. 2. If the joystick is in the initial position but still emits seven alarms, it indicates a malfunction in the controller circuit.
Eight short beeps	System failure in the controller.	Please consult Table 1-2. Please consult your distributor for maintenance.
Nine short beeps	Controller overload.	Controller overload; wheelchairs have excessively high inclinations or heavy loads when climbing slopes. After turning off the controller and restarting it, the failure should be resolved.
Ten short beeps	Battery failure.	3. The battery power is very low. Please, charge the battery. 4. Please check if the plug connected to the battery is normal and if the battery output voltage is normal.
Eleven short beeps	Controller overheating.	The controller is overheated. Turn off the controller and allow it to cool for a period of time before using it again.

6. Controller

6.7 Controller (Joystick)


Power Button

To turn on the wheelchair, press the power button. To turn it off, press the button again. When the wheelchair is on, the LCD screen will light up.

Sleep Mode

If the joystick is not used for more than 20 minutes, the power will automatically turn off, and the system will enter sleep mode.

The system can be turned on again by pressing the power button.

Speed Adjustment

The driving speed of the wheelchair can be adjusted according to the user's preference and circumstances.

Use the decrease and increase buttons to adjust the speed as needed.

Horn Button

Press the horn button to sound the horn.

Mute Button

Press the mute button to enable/disable the sound.

SOS Button

Press the SOS button to emit an emergency sound.

Battery Power Indicator

The battery indicator on the LCD screen is active when the power is on and shows the remaining battery power.

The indicator levels are as follows:

- 5 Bars: 80-100% remaining battery power
- 4 Bars: 60-80% remaining battery power
- 3 Bars: 40-60% remaining battery power


- 2 Bars: 20-40% remaining battery power
- 1 Bar: Less than 20% remaining battery power

When only one bar remains, it indicates low power, meaning the batteries need to be recharged.

Always ensure that the batteries are fully charged before taking long trips where recharging may not be available.

Do not operate the unit if only the single red bar is displayed, as this indicates an immediate need for recharging.

If the battery indicator blinks, it means less than 5% of the battery power remains.


Using the Joystick

- To control the speed and direction of the wheelchair, use the joystick.
- The distance you push the joystick determines the movement speed.
- Always ensure the joystick is in the central position when turning the unit on or off. If the joystick is not centered when turning the unit on or off, an audible alert will sound. To mute the alert, return the joystick to the center.

WARNING! In case of an emergency, you can also turn off the power directly by pressing the power button.

6.8 Using the Remote Control

Battery Insertion and Remote Control Connection

- Insert AAA batteries into the remote control.
- Turn on the remote control by pressing the "power button" for 3 seconds.
- You will know the remote control is on when the green lights on the power button start flashing.

Power Button

To turn on the remote control, press the power button.

To turn it off, press the button again.

Note: The remote control does not have the ability to turn your wheelchair on or off. It is designed to provide the ability to control your wheelchair remotely through the joystick control. To use your remote control, the main control must be turned on.

Speed Adjustment

The wheelchair's driving speed can be adjusted remotely using the remote control to match the user's preference and circumstances. Simply use the decrease and increase buttons to adjust the speed as needed.

Horn Button

Press the horn button to sound the horn.

Using the Joystick (Wheelchair Direction)

To control the speed and direction of the wheelchair, use the joystick.

The distance you push the joystick determines the movement speed.

Automatic Recline Button

- Press and hold the backrest button (B) until the button lights up.
- Press and hold the acceleration or deceleration button to recline or straighten the wheelchair.
- Once the reclining process is complete, briefly press button B to return to driving mode.

Auto-Unfold (Open the Wheelchair) and Auto-Fold (Close the Wheelchair) Buttons

- Press and hold the fold button for 3 to 5 seconds; the folding indicator will light up.
- Press and hold the acceleration or deceleration button to fold or unfold the wheelchair.
- Once the folding or unfolding process is complete, briefly press the fold button to return to driving mode.

WARNING! In an emergency, you can also turn off the power directly by pressing the power button.

The function buttons on the remote control are designed for specific versions, and not all wheelchair versions can use them.

6.9 Remote Control Synchronization

In general, the remote control has been successfully paired with the controller in the same package before shipping, so there is no need to pair it again.

If the remote control cannot connect to the controller or when purchasing a new remote control, follow these methods:

Synchronizing Your Remote Control

Battery Insertion and Turning on the Remote Control

- Insert AAA batteries into the remote control.
- Turn on the remote control by pressing the "power button" for 3 seconds.
- You will know the remote control is on when the green lights on the power button start flashing.

Starting Synchronization Mode

- To pair your remote control with the wheelchair, press the "deceleration button" on the remote control for 3 seconds.
- The green light on the power button will turn red, indicating that the pairing search is activated.

Synchronizing with the Joystick Controller

- Turn on the joystick controller.
- Press and hold the deceleration button on the controller for 3 seconds. The pairing red light will flash.
- Press and release the horn button on the controller. The pairing red light should stop flashing and remain lit, indicating a successful synchronization.

- If the remote control has not been successfully paired, repeat the pairing process.

Completing the Pairing Process

- After pairing, restart the controller by pressing the power button.
- Press the power button on the remote control to exit pairing mode. The green power light on the remote control will flash.
- The pairing indicator light on the controller will remain lit, indicating that the remote control is now ready to operate the wheelchair.

7. Maintenance

Wheelchair maintenance includes cleaning the wheelchair, checking the wheels and battery, and charging the battery. For additional maintenance, contact your distributor. We suggest checking the wheelchair every six months and performing an annual inspection.

7.7 Battery Cleaning and Maintenance

Wheelchair Cleaning

Please clean the wheelchair regularly.

Wipe down parts frequently touched by the user (such as the seat, armrests, and controller) with a clean, slightly damp cloth. Do not use organic solvents for cleaning.

If the user is a patient, the wheelchair should be cleaned once a week.

If the wheelchair is used by a person with an infection, it should be cleaned and disinfected with disinfectants.

- **Wheels**

Regularly check tire pressure and wear condition. When the tire tread depth is reduced to 1 mm, replace it with a new tire.

- **Battery**

Ensure the battery is fully charged frequently. To extend the battery life, we strongly recommend that users do not let the battery fully discharge before recharging.

7.8 Wheelchair Malfunction and Inspection

If malfunctions occur during wheelchair operation, turn off the power before inspecting.

Symptom: Complete power loss, and all LED lights on the controller panel are off.

Inspection Steps:

Step 1: Check if the controller plug is loose.

Step 2: Check if the connection between the controller plug and the battery box is loose. Please reinsert the plug connector (hold the plug when removing it. Do not pull the cable to avoid unnecessary damage to the wiring).

If, after the above checks, the wheelchair still does not regain power, or if the user has any doubts about the inspection steps, please contact your distributor.

The controller has a diagnostic system to monitor the controller and motor. Any malfunctions

in these parts will be indicated by the controller. For more details, refer to the Audio Signal Indication Chapter.

7.9 Maintenance Inspection

The following is a checklist. The electric wheelchair should be inspected according to our recommendation.

When entering or exiting the wheelchair, some automatic checks are performed.

For more attention, we have specially listed these self-check items in Area A.

12. Warranty

- **Warranty Coverage**

We have carefully designed this wheelchair for you.

If there are defective materials or manufacturing errors, we provide free repairs and lifetime maintenance under the terms and conditions set forth in the warranty cards.

- **The warranty does not cover the following conditions:**

Subjective evaluations without functional problems.

Wear and aging (coating surface and electroplating, natural resin discoloration, etc.).

- **Maintenance is not applicable in the following cases:**

Failure to regularly perform the checks specified by us.

Improper or incorrect maintenance.

Operation different from that indicated in the manual or overloading.

Unauthorized modifications.

External factors, such as soot, pharmaceuticals, bird droppings, acid rain, flying stones, metal dust, etc.

Natural disasters, such as typhoons, floods, fires, earthquakes, etc.

- **The following costs are not covered:**

Replacement of consumable materials, such as tires, fuses, plastic parts, glass parts, lubricants, etc.

Inspection, adjustment, oil addition, cleaning fees, etc.

Costs for regular checks specified by our factory.

Unauthorized modifications.

Maintenance costs outside our authorized service centers.

- **Customer Precautions**

To ensure the warranty remains valid, customers are required to follow these precautions:

Test and properly drive the electric wheelchair according to the manual instructions.

Perform daily inspections.

Implement verification recommendations according to our suggestion.

- **Warranty Acceptance**

If warranty service is required, bring the electric wheelchair and present the after-sales service card at our service center for warranty service. If the user cannot provide the after-sales service card, the service will not be provided.

- **Environmental Protection**

To protect the environment, any damaged or discarded parts of the electric wheelchair should be returned to our factory or delivered to the State Department for disposal. Do not dispose of them improperly.

- **Warranty Validity**

The after-sales service card is valid from the date stamped.

This product warranty is valid for 1 year under normal operating conditions.

Spare parts (such as battery pack, seat cushion, tires, side panel, armrest cushion, etc.) are not covered under the warranty.

During the warranty period, maintenance is free.

After the warranty period, maintenance will be charged; exact costs should be consulted at the service center.

13. Manufacturer's Guidelines and Declaration

Electromagnetic Emissions:

Manufacturer's Guidelines and Declaration - Electromagnetic Immunity			
Powered Wheelchairs must be used in the specified electromagnetic environment described below, and the buyer or user of the powered wheelchair must ensure that the following conditions are met:			
Anti-Interference Measurement	IEC60601 Test Level	Compatibility Level	Electromagnetic Environment Conditions
Electrostatic Discharge (ESD) ISO7176 IEC60601	±6 kV Contact Discharge ±8 kV Air Discharge	±6 kV Contact Discharge ±8 kV Air Discharge	The floor should be made of wood, concrete, or ceramic, and if covered with synthetic materials, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst ISO7176 IEC60601	±1 kV Voltage in Differential Mode ±2 kV Voltage in Common Mode	±2 kV To the Power Supply Cable	The power supply in hospitals or commercial environments must be of standard quality.
Surge ISO7176 IEC60601	±2 kV To the Power Supply Cable	±2 kV To the Power Supply Cable	The power supply in hospitals or commercial environments must be of standard quality.
Voltage dips, short interruptions, and voltage variations at power input	<5% UT, Duration of 0.5 cycles (at UT, >95% drop) 70% UT, Duration of 25	<5% UT, Duration of 0.5 cycles (at UT, >95% drop)	The power supply in hospitals or commercial environments must be of

	cycles (at UT, 30% drop) <5% UT, Duration of 5 seconds (at UT, >95% drop)	40% UT, Duration of 1 cycle (at UT, 100% drop) 70% UT, Duration of 25 cycles (at UT, 30% drop) <5% UT, Duration of 5 seconds (at UT, >95% drop)	standard quality. If a powered wheelchair user needs to operate the wheelchair during a power outage, battery power is recommended.
Power Frequency Magnetic Field (50/60Hz)	30 A/m	30 A/m	The power frequency magnetic field should have the typical horizontal characteristics of a power frequency magnetic field in a commercial or hospital environment.
Note: UT refers to the AC mains voltage before the test voltage is applied.			

Motorized Wheelchairs Must Be Used in the Electromagnetic Environment Specified Below The buyer or user of the motorized wheelchair must ensure that the following conditions are met:			
Anti-Interference Measurement	Test Level IEC60601	Electromagnetic Compatibility Level	Environmental Conditions
RF Conducted Disturbances ISO7176 IEC60601	3 Vrms From 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communication equipment should not be used near any part of the motorized wheelchair, including cables. Respect the recommended isolation distance. The distance should be calculated according to the corresponding formula for the transmitter frequency.
RF Radiation (Charger) ISO7176 IEC60601	3 V/m From 80 MHz to 1.0 GHz	3 V/m	Recommended Isolation Distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ (80 MHz to 800 MHz) $d = 2.3 \sqrt{P}$ (800 MHz to 1.0 GHz) Recommended Isolation Distance $d = 0.2 \sqrt{P}$ (26 MHz to 800 MHz) $d = 0.4 \sqrt{P}$ (800 MHz to 2.5 GHz)
RF Radiation (Wheelchair) ISO7176 IEC60601	20 V/m From 26 MHz to 25 GHz	20 V/m	Where P is the maximum nominal output power of the transmitter provided by the manufacturer in watts (W), and d is the recommended isolation distance in meters (m). The intensity of the fixed RF transmitter field is determined by electromagnetic field measurement A, in each frequency band. B should be lower than the compatibility level. Interference may occur near equipment

			containing the following symbols: 
OBSERVATION 1: At a frequency of 80 MHz and 800 MHz, the higher frequency range applies.			
OBSERVATION 2: These guidelines may not be suitable for all cases, for example, when electromagnetic propagation is affected by the absorption and reflection of buildings, objects, and human bodies.			

The motorized wheelchairs must be used in the specified electromagnetic environment below, and the buyer or user of the motorized wheelchair must ensure compliance with the following conditions:			
Anti-Interference Measurement	Test Level IEC60601	Electromagnetic Compatibility Level	Environmental Conditions
RF Conducted Disturbances (ISO7176, IEC60601)	3 Vrms From 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communication equipment should not be used near any part of the motorized wheelchair, including cables. Respect the recommended isolation distance. The distance should be calculated according to the corresponding formula for the transmitter frequency.
RF Radiation (Charger) ISO7176, IEC60601	3 V/m From 80 MHz to 1.0 GHz	3 V/m	Recommended Isolation Distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ (80 MHz to 800 MHz)}$ $d = 2.3 \sqrt{P} \text{ (800 MHz to 1.0 GHz)}$ Recommended Isolation Distance $d = 0.2 \sqrt{P} \text{ (26 MHz to 800 MHz)}$ $d = 0.4 \sqrt{P} \text{ (800 MHz to 2.5 GHz)}$
RF Radiation (Wheelchair) ISO7176, IEC60601	20 V/m From 26 MHz to 25 GHz	20 V/m	Where P is the maximum nominal output power of the transmitter provided by the manufacturer in watts (W), and d is the recommended isolation distance in meters (m). The intensity of the fixed RF transmitter field is determined by electromagnetic field measurement A, in each frequency band. B should be lower than the compatibility level. Interference may occur near equipment containing the following symbols: 
OBSERVATION 1: At a frequency of 80 MHz and 800 MHz, the higher frequency range applies.			
OBSERVATION 2: These guidelines may not be suitable for all cases, for example, when electromagnetic propagation is affected by the absorption and reflection of buildings, objects, and human bodies.			

Post-Sales Service Card

Cardholder:	Contact Phone:	Model:	
Add:		Manufacturing Date:	
Weight:	Height:	Joystick Position:	Serial Number:
Quality Warranty: Your purchased product includes the following warranty services: 1. The motor, controller, and structure have a one-year warranty. During the warranty period, maintenance is free of charge. 2. We provide door-to-door maintenance; service charges must be covered according to the distance.			
Service Record	Date	Description	

Note:

1. If you have changed your contact information, please inform us in time.
2. Replaced parts are the property of our company.
3. Consumable parts such as batteries, tires, armrests, etc., are not covered by the warranty.

CONTACT INFORMATION:

Address: Rua Brusque, 1078, Ground Floor Box 02,

Centro, Postal Code: 88.302-001, Itajaí-SC

CNPJ: 45.784.439/0001-90

Technical Responsible: Priscila de Gouveia Zamboni Lins CRF-SC 20958

Versão 03

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Registro ANVISA: 82692690005

Versión 03

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Version 03

Date: December 2024